














Equipment	Recommended Replacement Schedule	
CPAP (Single Pressure Device)	1 per 5 years	
Auto CPAP (Adjusting Pressure)	1 per 5 years	Same as above with different internal mechanisms and options
BiLevel S – (Bi-Level Pressure)	1 per 5 years	
Auto Bi-Level (Adjusting Bi-Level Pressure)	1 per 5 years	Same as above with different internal mechanisms and options
Full Face Mask	1 per 3 months	
Full Face Mask Cushion Replacement	1 per month	

<p>Nasal Mask</p>	<p>1 per 3 months</p>	
<p>Nasal Mask Cushion Replacement</p>	<p>2 per months</p>	
<p>Nasal Pillow Mask/System</p>	<p>1 per 3 months</p>	
<p>Nasal Pillow Replacements</p>	<p>2 per month</p>	
<p>Combination Oral/Nasal Mask</p>	<p>1 per 3 months</p>	

<p>Oral Interface</p>	<p>1 per 3 months</p>	
<p>Headgear</p>	<p>1 per 6 months</p>	
<p>Chinstrap</p>	<p>1 per 6 months</p>	
<p>CPAP Tubing</p>	<p>1 per 3 months</p>	
<p>Climateline Tubing for ResMed S9</p>	<p>1 per 3 months</p>	
<p>Disposable Filter</p>	<p>2 per month</p>	

Non-Disposable Filters	1 per 6 months	
Comfort Care Pad*	1 per month	
Heated Humidifier	1 per 3 years	
Replacement Humidifier Chamber	1 per 6 months	